




SOLACE[®] PREVENTION

SERIES THERAPEUTIC SUPPORT SURFACES




The **Invacare America Solace[®] Prevention Series Mattresses** are unique foam mattress replacements that are designed to help minimize the risk of pressure injuries by distributing the weight as uniformly as possible to allow for pressure redistribution. The unique and patented designs in the Prevention Series combine comfort and technology for resident and caregiver confidence.

Invacare America Solace® Prevention 1080 Mattress

The convoluted “fortress” cut design offers more surface area for resident support and relief. This high quality foam includes a soft head and heel section. Available in 80" standard length or custom 84" (SPS1084) and optional raised side rail design (SPS1080RSR). The three zones help relieve pressure in the important susceptible body zones.

Features

- 70 Denier nylon top
- 11 oz. tri-laminate vinyl bottom
- Fluid resistant
- Raised side rail design available
- Pressure reduction
- Compliant with federal fire regulations*




Invacare America Solace® Prevention 2080 Mattress

Added to the design are built-in side supports. The head and heel soft sections add comfort and help with pressure redistribution, while the side supports provide stability and durability. Also available in bariatric sizes; SPS2080B42, SPS2080B48 (450 lb. weight capacity on bariatric models).

Features

- Built-in side supports
- 70 Denier nylon top cover
- 11 oz. tri-laminate vinyl bottom
- Soft head and heel section
- Fortress cut convoluted foam
- Pressure reduction



1. Test shows subject in supine position demonstrating contact measurement between key pressure points using BodiTrak mapping technology on the SPS1084 support surface. Subject: Male 5'11", 185 lb.

Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
SPS1080	80" L 36" W 6" H	350lb.	16 CFR1632 16 CFR1633	5 Years (Mattress) 1 Year (Cover)
SPS1084	84" L 36" W 6" H	350lb.	16 CFR1632 16 CFR1633	5 Years (Mattress) 1 Year (Cover)
SPS1080RSR	80" L 36" W 6" H 3" Raised Side Rails	350lb.	16 CFR1632 16 CFR1633	5 Years (Mattress) 1 Year (Cover)

Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
SPS2080B42/ SPS2080B48	80" L 42"/ 48" W 6" H	400lb.	16 CFR1633 16 CFR1633	5 Years (Mattress) 1 Year (Cover)

*Federal flammability standards: 16CFR1632, 16CFR1633, CA Technical Bulletin 129