

Features




Technical data

For more information about this product, including the user manual, please visit our website: www.invacareamerica.com

	Seat width	Overall width	Folded width	Seat depth	Overall depth (without riggings)	Overall height
Transport Classic	19 inches	22 inches	9 inches	16 inches	18 inches	39.5 inches

	Seat to floor height	Back height	Rear wheel size	Weight	Shipping weight	Weight limit
Transport Classic	19.5 inches	18.5 inches	8 inches (All Wheels)	28.88 lbs	31.63 lbs	250 lbs

Item No.	Description
EGTC19STFSBK	Transport Classic 19", Fixed Arms, Swing Away Footrests – Black
EGTC19STFSBU	Transport Classic 19", Fixed Arms, Swing Away Footrests – Blue



HCPCS Code: Pending

EasyGlide™ Transport Classic



Invacare America
One Invacare Way
Elyria, Ohio 44035-4190
(800) 333-6900

©2026 Invacare America. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare® America or its subsidiaries unless otherwise noted. Specifications are subject to change without notification. Form No. 26-002, Job 260105



Durable performance
for dependable caregiver-assisted transport



Introducing the new EasyGlide™ Transport Classic

A dependable, caregiver-assisted transport chair designed for safe mobility.

The EasyGlide™ Transport Classic comes fully assembled and is designed for caregiver-assisted mobility in everyday settings.

Featuring fixed arms and a durable steel frame, it delivers dependable transport support across care environments.

Engineered for stability and straightforward handling, the Classic includes standard swing-away footrests and a folding backrest to simplify transfers and storage.




Reliable Performance

A reinforced frame supports up to 250 lbs and is suited for routine transport use. Integrated wheel locks help secure the chair during transfers and stationary positioning.



Supportive Comfort

Seat and back padding provide cushioned support during transport, while padded armrests enhance comfort and stability. The seating surface is designed to support upright positioning for short-term mobility needs.



Foldable & Compact

The backrest folds down to reduce overall size when not in use. Measuring just 9 inches wide when folded and 29 inches high from floor to top of chair, its compact profile and manageable weight allow for easier lifting and vehicle placement.

Ergonomic push handles

Storage pocket

Seatbelt for added safety

Push-to-lock wheel brakes

Integrated curb assistants

Fold-down backrest

For compact storage

Padded seat, back and arms

Steel frame

Swing away foot rests

High-quality 8" TPU casters

*Product specifications may vary from those displayed.

3 Year Limited Warranty