

Features



Comfortable and padded



Compact and travel ready



Fold down backrest



Removable, swing away arms

Technical data

For more information about this product, Including the user manual, please visit our website: www.invacareamerica.com

EasyGlide™ Wheelchair	Seat width	Overall width	Folded width	Seat depth	Overall depth (without riggings)	Overall height
	16 inches 18 inches 20 inches	29.5 inches	14.5 inches	17.72 inches	Not Applicable	36 inches
EasyGlide™ Wheelchair	Seat to floor height	Back height	Rear wheel size	Weight	Shipping weight	Weight limit
	19.5 inches	18.5 inches	24 inches	41.33 lbs	47.83 lbs	300 lbs

Tested to meet ISO 7176 and RESNA standards for strength, durability and safety

Item No.	Description
EGWC16STRSBU	16", Removeable Desk Arms, Swing Away Footrests – Blue (Steel)
EGWC16STRSBK	16", Removeable Desk Arms, Swing Away Footrests – Black (Steel)
EGWC18STRSBK	18", Removeable Desk Arms, Swing Away Footrests – Black (Steel)
EGWC18STRSBU	18", Removeable Desk Arms, Swing Away Footrests – Blue (Steel)
EGWC20STRSBK	20", Removeable Desk Arms, Swing Away Footrests – Black (Steel)
EGWC20STRSBU	20", Removeable Desk Arms, Swing Away Footrests – Blue (Steel)
EGWC16STREBK	16", Removeable Desk Arms, Elevated Legrests – Black (Steel)
EGWC16STREBU	16", Removeable Desk Arms, Elevated Legrests – Blue (Steel)
EGWC18STREBK	18", Removeable Desk Arms, Elevated Legrests – Black (Steel)
EGWC18STREBU	18", Removeable Desk Arms, Elevated Legrests – Blue (Steel)
EGWC20STREBK	20", Removeable Desk Arms, Elevated Legrests – Black (Steel)
EGWC20STREBU	20", Removeable Desk Arms, Elevated Legrests – Blue (Steel)

Invacare America
One Invacare Way
Elyria, Ohio 44035-4190
(800) 333-6900

©2026 Invacare America. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare® America or its subsidiaries unless otherwise noted. Specifications are subject to change without notification. Form No. 26-003, Job. 260106

Black

Blue

Frame Color Frame Color

HCPCS Code: Pending



EasyGlide™ Wheelchair



Durable performance
for independent everyday use



Introducing the new EasyGlide™ Wheelchair

A reliable and robust manual wheelchair designed to support everyday mobility.

Built with durable steel construction and tested for strength and safety, the EasyGlide™ Wheelchair delivers dependable performance wherever you go. Designed for ease of use and portability, the EasyGlide™ Wheelchair folds compactly and stores with ease.

Flip-back, removable desk arms and the option of swing-away or elevating leg rests facilitate easier transfers and flexibility for daily activities.

Shipped fully assembled, the EasyGlide™ Wheelchair is ready to use right out of the box.




Reliable Performance

Constructed with a sturdy steel frame and solid flat-free tires, the EasyGlide™ Wheelchair provides stability and durability for daily use. The push-to-lock braking system enhances safety so you feel secure and supported.



Configurable Support

Designed to adapt to individual needs, the EasyGlide™ Wheelchair offers multiple configuration options. Choose between swing-away or elevating leg rests, removable flip-back desk arms, and a range of seat widths to accommodate different users and care requirements.



Foldable & Compact

The fold-down backrest and streamlined frame help to ensure folding is simple and effortless. Compact when folded to just 12.5 inches wide, storage and transportation are made easy.

Ergonomic push handles

Padded seat and back

Storage pocket

Easy to reach push to lock wheel brakes

Integrated curb assistants

Fold-down backrest

For compact storage

Removable, swing away arms

Swing away or elevating leg rests available



*Product specifications may vary from those displayed.

5 Year Limited Warranty