

Features



Compact and travel ready



Comfortable and transfer friendly



Removable, swing away desk arms

Technical data

For more information about this product, including the user manual, please visit our website: www.invacareamerica.com

	Seat width	Overall width	Folded width	Seat depth	Overall depth (without riggings)	Overall height
						
Transport Flex	19 inches	22 inches	9 inches	16 inches	18 inches	39.5 inches

	Seat to floor height	Back height	Rear wheel size	Weight	Shipping weight	Weight limit
						
Transport Flex	19.5 inches	18.5 inches	8 inches (All Wheels)	Steel: 36.37 lbs Aluminum: 25.13 lbs	42.87 lbs	250 lbs

Item No.	Description
EGTC19STRSBK	Transport Flex 19", Removeable Desk Arms, Swing Away Footrests – Black (Steel)
EGTC19STRSBU	Transport Flex 19", Removeable Desk Arms, Swing Away Footrests – Blue (Steel)
EGTC19ALRSBK	Transport Flex 19", Removeable Desk Arms, Swing Away Footrests – Black (Aluminum)
EGTC19ALRSBU	Transport Flex 19", Removeable Desk Arms, Swing Away Footrests – Blue (Aluminum)

Black

Frame Color

Blue

Frame Color

HCPCS Code: Pending

Invacare America
One Invacare Way
Elyria, Ohio 44035-4190
(800) 333-6900

©2026 Invacare America. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare® America or its subsidiaries unless otherwise noted. Specifications are subject to change without notification. Form No. 26-002, Job 260105



EasyGlide™ Transport Flex



Durable performance
for dependable caregiver-assisted transport



Introducing the new EasyGlide™ Transport Flex

A dependable, caregiver-assisted transport chair designed for safe mobility.

The EasyGlide™ Transport Flex comes fully assembled and is designed for caregiver-assisted mobility in everyday settings. It features swing-away, removable desk arms that swing upward for transfers or detach entirely, along with a steel or lightweight aluminum frame option for adaptable transport support.

Built for portability and controlled handling, the Flex includes standard swing-away footrests and a folding backrest to streamline transfers and storage.




Reliable Performance

A reinforced frame supports up to 250 lbs and is suited for routine transport use. Integrated wheel locks help secure the chair during transfers and stationary positioning.



Supportive Comfort

Seat and back padding provide cushioned support during transport, while the flip-back, removable desk arms offer both comfort and flexibility. The seating surface supports upright positioning during short-term caregiver-assisted mobility.



Foldable & Compact

The backrest folds down to reduce overall size when not in use. Measuring just 9 inches wide when folded and 29 inches high from floor to top of chair, its compact profile and manageable weight allow for easier lifting and vehicle placement.

Ergonomic push handles

Storage pocket

Seatbelt for added safety

Push-to-lock wheel brakes

Integrated curb assistants

Fold-down backrest

For compact storage

Removable, swing away arms

Steel or lightweight aluminum frame available

Swing away foot rests

High-quality 8" TPU casters

*Product specifications may vary from those displayed.

3 Year Limited Warranty